

Bag Features -

- A non-prescriptive design which can be used by anyone delivering different types of care.
- Easily adjustable storage system enables the bag to carry a wide range of equipment.
- Ergonomic design with padded shoulder straps and back.
- The inside of the bag is lined with easy to clean fabric.
- Strong durable material used throughout the bag.
- Free standing design.
- Tamper proof and lockable zips with easy to use large zip pulls.
- Integrated carry handle.
- Identification panel.
- Mesh side pockets for hand gel or drinks bottles.
- Separate paper work and laptop storage.

BISTON®

CONTAINING AND
LIGHTENING THE LOAD



Who we are...

Biston produces British made ergonomic products which are designed to lighten the load and address issues of load handling. The brand's two current products the Biston Point of Care Bag and the Biston Wheelchair Bag were both developed by Sue Barton and Jean Bish who are experienced chartered physiotherapists.

Biston is part of Sutherland Health Ltd, a supplier of innovative health care brands and products to the UK National Health Service (NHS) and other UK customers.

Sutherland Health Ltd is ISO approved and MHRA registered and has quality at the core of its business.

Visit our website at:

www.sutherlandhealth.com
email us at: info@sutherlandhealth.com
Call us: +44(0)1635 874488

Write to us:

Sutherland Health Ltd,
Unit 1, Rivermead, Pipers Way,
Thatcham, Berkshire, RG19 4EP.UK

The Biston® Point of Care Bag

The complete equipment storage system




Designed for those delivering
care in the community

The Biston Point of Care Bag has been specially designed by chartered physiotherapists to effectively store a wide range of equipment, while spreading the weight evenly to provide a comfortable experience and therefore lightening the load and making it safer to carry equipment.

A large front storage section is a key feature of the Biston Point of Care Bag. This section includes adjustable dividers, which can be fastened to secure a wide range of equipment. The bag also features a detachable pouch which can be used to separately store small items including waste material and sharps.

Easy to clean material is used to line the inside of the Biston Point of Care bag and antibacterial fabric is used for the detachable examination mat which can be neatly stored in the bag when not in use. The bag has a back office section which can be used to store pens and paperwork, it is also large enough to fit laptops and electronic tablets.



Bag Dimensions:

Weight (inclusive of all inserts)	2 kgs
Height (top/bottom)	49 cm
Width (at widest point)	33 cm
Depth (overall without divider)	23.5 cm
Depth (office section)	5 cm